



Michael R. Pence  
Governor

William C. VanNess II, MD  
State Health Commissioner

**DATE:** September 16, 2014  
**TO:** All Local Health Departments  
Attn: Chief Food Inspection Officer  
**FROM:** *Laurie Kidwell*  
Laurie Kidwell, RRT Supervisor  
Food Protection Program  
**SUBJECT:** Energizer Personal Care – Playtex Manufacturing, Inc. - RECALL[Medical Device]

**\*\*\*Expanded Recall\*\*\***

**AFFECTED PRODUCT:** Playtex Nurser Deluxe Double Electric Breast Pump

**SUMMARY:** Unclassified Recall; In March 2014, Playtex announced a voluntary recall of certain AC/DC power adapters because the casing on some adapters may come loose and separate, resulting in a potential for electric shock. After further evaluation and investigation, Playtex has added one adapter production lot to the recall.

The affected adapters were manufactured from November 2012 through July 2013. The products can be identified by product serial number (P12315 -XXXX through P13205-XXXX). Alternatively, the product can be identified by adapter production code (1238 through 1324). The adapters were sold along with the Playtex Nurser Deluxe Double Electric Breast Pump. They were not sold separately.

The Playtex Nurser Deluxe Double Electric Breast Pump was sold at nationwide, specialty and online retailers.

**SUGGESTED ACTION:** For consumer inquiry only. Consumers who have purchased an affected product should contact Playtex for a replacement by calling 1-888-207-1429 from 8 a.m. to 6 p.m. ET Monday through Friday or online at [www.playtexproducts.com](http://www.playtexproducts.com). Consumers should immediately discontinue use of the adapter if it shows signs of separating.

\*\*\*\*\*

**Recall -- Firm Press Release**



2 North Meridian Street • Indianapolis, IN 46204  
317.233.1325 tdd 317.233.5577  
[www.statehealth.in.gov](http://www.statehealth.in.gov)

To promote and provide  
essential public health services.

FDA posts press releases and other notices of recalls and market withdrawals from the firms involved as a service to consumers, the media, and other interested parties. FDA does not endorse either the product or the company.

## **Playtex Adds Production Lot To Its March 2014 Voluntary Recall Of Certain AC/DC Power Adapters Used With the Playtex® Nurser Deluxe Double Electric Breast Pump**

**Contact:**

Consumer:

888-207-1429

[www.playtexproducts.com](http://www.playtexproducts.com)

**Media:**

Catherine McCormack

212-819-4816

**FOR IMMEDIATE RELEASE** - September 9, 2014 - **Energizer Personal Care – Playtex Manufacturing, Inc.** announced today that, out of an abundance of caution, it is adding an adapter production lot to its voluntary nationwide recall of certain AC/DC power adapters that are used with the Playtex® Nurser Deluxe Double Electric Breast Pump. No injuries have been reported to date.

In March 2014, Playtex announced a voluntary recall of certain AC/DC power adapters because the casing on some adapters may come loose and separate, resulting in a potential for electric shock. After further evaluation and investigation, Playtex has added one adapter production lot to the recall.

The affected adapters were manufactured from November 2012 through July 2013. The products can be identified by product serial number (**P12315 -XXXX through P13205-XXXX**). Alternatively, the product can be identified by adapter production code (**1238 through 1324**). The adapters were sold along with the Playtex Nurser Deluxe Double Electric Breast Pump. They were not sold separately.

The Playtex Nurser Deluxe Double Electric Breast Pump was sold at nationwide, specialty and online retailers. Playtex is notifying its retail partners to return any remaining products with affected AC/DC adapters.

Consumers who have purchased an affected product should contact Playtex for a replacement by calling 1-888-207-1429 from 8 a.m. to 6 p.m. ET Monday through Friday or online at [www.playtexproducts.com](http://www.playtexproducts.com). Consumers should immediately discontinue use of the adapter if it shows signs of separating.

This recall is specific to certain AC/DC adapters sold with the Playtex Nurser Deluxe Double Electric Breast Pump, and is being conducted in conjunction with the U.S. Food and Drug Administration (FDA).

Adverse reactions or quality problems experienced with the use of this product may be reported to the FDA's MedWatch Adverse Event Reporting program either online, by regular mail or by fax.

- Complete and submit the report **Online**: [www.fda.gov/medwatch/report.htm](http://www.fda.gov/medwatch/report.htm)
- **Regular Mail or Fax**: Download form [www.fda.gov/MedWatch/getforms.htm](http://www.fda.gov/MedWatch/getforms.htm) or call 1-800-332-1088 to request a reporting form, then complete and return to the address on the pre-addressed form, or submit by fax to 1-800-FDA-0178

---

---

## Playtex Adds Production Lot To Its March 2014 Voluntary Recall Of Certain AC/DC Power Adapters Used With the Playtex® Nurser Deluxe Double Electric Breast Pump Photos



*Playtex*<sup>®</sup>

1250

I. T. E. POWER SUPPLY

MODEL: SY-12120


INPUT: 120V AC 60Hz 200mA

OUTPUT: 12V DC 1200mA


POLARITY: 



LISTED  
E221142  
1NW6



MADE IN CHINA



E239834  
14CC

**LISTED**

BREAST PUMP  
P/N: PYX-020  
ELEC. RATING: 12V DC 3200 mA  
MADE IN CHINA

S/N: R12206-759520



